



Fisher Hamilton®

Prestige™

Contemporary, yet timeless in its design, Prestige sets a new standard in steel laboratory casework.



Prestige combines soft, radiused fronts with concealed hinges, stylish pull options and a choice of 18 powder coat colors to create its look of distinction.

The strength and durability of Prestige's steel construction ensures that both the look and the function of your lab will endure.

Prestige... a new standard of form and function



Durable Steel Cabinet Construction

The Prestige steel cabinet is constructed of a one-piece wrap around cabinet shell with double corner posts in all four corners. The back of the cabinet is available with or without a service chase access panel. This durable construction will support up to 500 lbs. per lineal foot. Levelers are accessed through the toe base panel, which leaves the cupboard bottom perforation free. One piece, die-formed steel drawer bodies are designed for the rigors of the laboratory. Prestige is finished with a chemical and abrasion resistant urethane powder-coated finish for maximum durability.

Stylish Pulls



Broad, Comprehensive Product Offering

Prestige is available in floor mounted, suspended and mobile cabinets in a wide variety of door and drawer configurations and cabinet styles. A full range of work surfaces, fixtures and accessories are available to complete any project.

Self-closing, Concealed Hinges



18 Standard Colors

