

Mobile Equipment Cabinets

FROM FISHER HAMILTON

The Fisher Hamilton Mobile Equipment Cabinets are designed with a special high impact PVC edge around all doors, drawers, sides and countertops. The PVC banding absorbs impact and seals the edge so dirt, dust, water or chemicals cannot contaminate or damage the panel core material.



Mobile Equipment Cabinet

FEATURES

- Double doors that open 120 degrees for full cabinet access
- Six-inch high drawer with full-extension suspension
- Front pivotal, locking casters
- Height-adjustable inner shelf for storage
- Chemical resistant work surface
- Grommet cutouts for electrical connections

Product Number	Height	Depth	Width
955P275 D	37"	29"	30"

Double-Wide Mobile Equipment Cabinet

FEATURES

- Left-hand door that opens 120 degrees for full cabinet access
- Two six-inch high drawers with full-extension suspension
- Front pivotal, locking casters
- Pull-out shelf for storage of equipment
- Chemical resistant work surface
- Grommet cutouts for electrical connections



Product Number	Height	Depth	Width
955P276 D	37"	30"	48"

Mobile Computer Carts

FROM FISHER HAMILTON



Mobile Computer Cart

FEATURES

- Double doors open 90 degrees, and then slide into the cabinet interior
- Pull-out keyboard drawer retracts into the cabinet
- Front pivotal, locking casters
- Three pull-out shelves for easy access to printer, server and paper supplies
- Chemical resistant work surface
- Grommet cutouts for electrical connections

Product Number	Height	Depth	Width
955P277 D	37"	29"	30"

Specifications

Substrate: High density 3/4" M-3 particle board with 1" plywood base for structural integrity.

Full Sub-Top: Full sub-top on all cabinets for strength and rigidity.

Worksurface: Worktops are manufactured from 1" particle board with horizontal grade laminates. All edges are banded with a high impact 3MM PVC edge banding.

Drawer Bodies: All drawers are manufactured with lock shoulder joinery. Exposed surfaces are melamine clad.

Drawers: Drawers have easy-to-open, 150 pound full extension 3 piece drawer suspension that will provide flawless performance over years of routine use.

Hinges: Five-knuckle institutional hinges are beveled at the top and bottom for easy cleaning.

Pulls: Recessed pulls are manufactured from injection molded plastic with a built in label holder.

Adjustability: Shelf and pull-out shelves are adjustable on 1.25" increments.

Locks: All doors and drawers are equipped with 5 disc tumbler locks.

Color Selection

All standard colors are WilsonArt selections in GP28 vertical surface grade laminate. Other WilsonArt selections, as well as those of other manufacturers, are available as non-standards.

Code	Color	Code	Color
1573-60-335	Frosty White (9Y)	D97-60-335	Haze (9S)
1570-60-335	White (LD)	D25-60-335	Atlantis (6F)
1572-60-335	Antique White (9P)	1503-60-335	Putty (9U)
1500N-60-335	Grey (9T)	D327-60-335	Pepperdust (6W)
D92-60-335	Dove Grey (9R)	1531-60-335	Light Beige (9X)
D381-60-335	Fashion Grey (7Y)	1511-60-335	Mandarin Red (9V)
D315-60-335	Platinum (6Q)	D414-60-335	Shogun (8N)
D30-60-335	Natural Almond (6I)	4622-60-335	Grey Nebula (LE)
D417-60-335	Lapis Blue (8Q)	4621-60-335	White Nebula (9Q)
D421-60-335	Midori (8S)	4628-60-335	Sand Nebula (LK)
D369-60-335	Burgundy (7Q)	4629-60-335	Blush Nebula (LL)
D455-60-335	Montpellier (LR)		
D72-60-335	Featherstone (8Y)		
1595-60-335	Black (90)		
D432-60-335	Cashmere (LQ)		
D378-60-335	Clear Teal (7W)		



Fisher Hamilton

Fisher Hamilton L.L.C.
1316 18th Street
Two Rivers, WI 54241
Tel: 920-793-1121
Fax: 920-793-3084
www.fisherhamilton.com

AL-1418-1 APR01
Copyright 2001 Fisher Hamilton L.L.C.