


SafeStor™ solvent storage cabinet features:

- Modular design
- Durable steel construction
- Heavy-duty load-bearing capacity
- Aesthetically attractive



SafeStor™ Solvent Storage Cabinets

Single Door Storage Cabinets




- Right hinged manual or self-closing door
- One adjustable shelf and spill tray at bottom
- Vent plug on left and right sides
- 17" (431mm) wide x 22" (559mm) high x 17" (432mm) deep
- Meets NFPA-30

Product No.	Type of Door	Capacity
950S8740YL	Manual	4 Gal.
950S8750YL	Self-Closing	4 Gal.



- Right hinged manual or self-closing door
- One adjustable shelf and spill tray at bottom
- Vent plug on left and right sides
- 23-1/4" (590mm) wide x 35" (889mm) high x 18" (457mm) deep


Product No.	Type of Door	Capacity
950S8760YL	Manual	12 Gal.
950S8770YL	Self-Closing	12 Gal.



- Right hinged manual closing door
- One adjustable shelf and spill tray at bottom
- Vent plug on left and right sides
- 23-1/4" (590mm) wide x 44" (1118mm) high x 18" (457mm) deep

Product No.	Type of Door	Capacity
950S8780YL	Manual	15 Gal.


Double Door Storage Cabinets – 36" Wide



- Double manual closing
- One adjustable shelf and spill tray at bottom
- Vent plug on left and right sides
- 36" (914mm) wide x 35" (889mm) high x 24" (610mm) deep


Product No.	Type of Door	Capacity
950S8670YL	Manual	30 Gal.

Double Door Storage Cabinets – 36" Wide (Cont.)



- Double manual or self-closing doors
- One adjustable shelf and spill tray at bottom
- Two vent plugs on rear of cabinet
- 36" (914mm) wide x 35-1/4" (897mm) high x 21-5/8" (549mm) deep


Product No.	Type of Door	Capacity
950S8800YL	Manual	22 Gal.
950S8810YL	Self-Closing	22 Gal.



- Double manual or self-closing doors
- Two adjustable shelves and spill tray at bottom
- Vent plug on left and right sides
- 36" (914mm) wide x 65" (1651mm) high x 34" (864mm) deep


Product No.	Type of Door	Capacity
950S8700YL	Manual	60 Gal.
950S8710YL	Self-Closing	60 Gal.

Double Door Storage Cabinets – 42" Wide



- Double manual closing doors
- Spill tray at bottom
- Vent plug on left and right sides
- 42" (1067mm) wide x 18-1/4" (464mm) high x 18" (457mm) deep


Product No.	Type of Door	Capacity
950S8720YL	Manual	12 Gal.



- Double manual closing doors
- One adjustable shelf and spill tray at bottom
- Vent plug on left and right sides
- 42" (1067mm) wide x 24-1/4" (617mm) high x 18" (457mm) deep


Product No.	Type of Door	Capacity
950S8730YL	Manual	17 Gal.

Double Door Storage Cabinets – 42" Wide (Cont.)




- Double manual or self-closing doors
- One adjustable shelf and spill tray at bottom
- Vent plug on left and right sides
- 42" (1067mm) wide x 44" (1118mm) high x 18" (457mm) deep

Product No.	Type of Door	Capacity
950S8650YL	Manual	30 Gal.
950S8660YL	Self-Closing	30 Gal.



- Double manual or self-closing doors
- Two adjustable shelves and spill tray at bottom
- Vent plug on left and right sides
- 42" (1067mm) wide x 65" (1651mm) high x 18" (457mm) deep

Product No.	Type of Door	Capacity
950S8680YL	Manual	45 Gal.
950S8690YL	Self-Closing	45 Gal.



- Double manual closing doors
- Five adjustable shelves and spill tray at bottom
- Vent plug on left and right sides
- 42" (1067mm) wide x 44" (1118mm) high x 12" (305mm) deep

Product No.	Type of Door	Capacity
950S8790YL	Manual	20 Gal.

General Information

- All of the SafeStor cabinets, except 950S8740 and 950S8750, have been submitted, tested and listed by FM (Factory Mutual) 6050.
- Dimensions are with levelers retracted.
- All shelves adjustable in nominal 1/2" (12mm) increments.
- Optional vents to fume hood piping assembly available when required.
- All solvent cabinets include locked doors.



Typical Detail – Spill Tray

Dimensions are nominal and illustrations and specifications are based on the latest product information available at the time of publication. The right is reserved to make changes at any time without notice.



Fisher Hamilton L.L.C.
1316 18th Street
Two Rivers, WI 54241
Phone: 920-793-1121
www.fisherhamilton.com